

ELMO

MG-600

MAGNISCOPE VIEWER FOR ELMO PRESENTERS



The Magniscope Viewer MG-600, used with ELMO EV/HV series Presenters, enlarges glass slides/medical slides and other miniature transparent objects placed on its object stage, and enables the enlarged image to be viewed on any monitor or through any video/data projector.

It is an extremely useful for scientific and educational communities.

The enlarged scientific view is a fun and great impression to all the students.



MG-600 on ELMO EV-4400AF Presenter


ELMO Magniscope Viewer MG-600

(Product Code # 9448)

Applicable Presenter	Magnifications (when viewed on 20" monitor)
EV-4400AF/4450AF	Approx. 200 times
EV-2000AF/2500AF (*)	Approx. 170 times
EV-400AF/450AF	Approx. 125 times
HV-5000XG/5500XG	Approx. 130 times
HV-3000XG/3500XG (*)	Approx. 150 times

Note : The optional light viewer panel is required for the use with the (*) marked Presenters.

PROFILE



ELMO ELMO CO., LTD.

6-14, Meizen-cho, Mizuho-ku, Nagoya, 467-8567, Japan E-Mail:foreign-div@elmo.co.jp

ELMO Mfg. Corp.

1478 Old Country Road,
Plainview, NY 11803-6034, U.S.A.
Tel. 516-501-1400 Fax. 516-501-0429
E-mail:elmo@elmousa.com
Web: <http://www.elmousa.com>

ELMO Canada Mfg. Corp.

44 West Drive, Brampton,
Ontario, L6T 3T6, Canada
Tel. 905-453-7680 Fax. 905-453-2391
E-mail:info@elmocanada.com
Web: <http://www.elmocanada.com>

ELMO (Europe) G.m.b.H.

Neanderstr 18,
40233 Düsseldorf, Germany
Tel. 0211-376051-53 Fax. 0211-376630
E-mail:elmo@europe.AOL.com
Web: <http://www.elmo.de>